

Tribin BA

Alkaline circulation cleaner with active chlorine

Properties

- Dissolves grease, protein and other organic soiling
- foam-free
- Especially for cleaning blood systems
- with water hardness stabilizers

Application

Tribin BA is an alkaline, chlorine-containing CIP and circulation cleaner for special use as a cleaning agent for blood-carrying systems in slaughterhouses. Tribin BA reliably cleans all areas of blood systems such as sticks, hoses and pipes, collection a

Technical data

Density (20°C)	pH-value
1,20	12 - 13 ; 1 %

Titration

Prepare 10 ml application solution and destroy the chlorine by adding a spatula tip of sodium thiosulphate. Titrate with 0.1 N hydrochloric acid against the indicator phenolphthalein.

Consumption (ml) x 0.27 = concentration in %

Notes

Store the product only in the original container and frost-free.
Direct sunlight should be avoided.
After removing partial quantities, the container should be resealed.

For commercial use only. This leaflet is for non-binding information only. The information is based on our current knowledge and experience. In any case, the user is obliged to carry out his own tests and trials to check the suitability of the products for his intended processes and purposes. The information in this leaflet does not constitute a guarantee for the quality and durability of the goods to be supplied by us. We reserve the right to make technical changes within the scope of what is reasonable. The current version of the corresponding EU safety data sheet must also be observed.