

# Mosal T

## Disinfectant concentrate

### Properties

- tested hygienic agent for water separator care in dry cleaning machines
- universally applicable liquid active ingredient concentrate with a pronounced broad-spectrum effect
- reliably combats mold, yeast and bacteria
- free of chlorine and formaldehyde

### Application

Water separator maintenance

After each water discharge, approx. 10 ml of Mosal T and approx. 100 ml of fresh water are added to the water separator. The previously emptied water separator is wiped out weekly with a cloth soaked in Mosal T.

### Technical data

Density (20°C)	pH-value
1,00 kg / l	6,0 - 8,0

### Biocide

Use biocide products with care. Always read the label and product information before use.

Notified according to national legislation for biocides under No. N-59810.

100 g product contains 2.45 g alkyl(C12-16)dimethylbenzylammonium chloride

Maximum usability from manufacture: 18 months in the closed original container.

### Notes

Store the product only in its original container and keep it frost-free.

After removing partial quantities, the container should be closed again.

For commercial use only. This leaflet is for non-binding information only. The information is based on our current knowledge and experience. In any case, the user is obliged to carry out his own tests and trials to check the suitability of the products for his intended processes and purposes. The information in this leaflet does not constitute a guarantee for the quality and durability of the goods to be supplied by us. We reserve the right to make technical changes within the scope of what is reasonable. The current version of the corresponding EU safety data sheet must also be observed.